

For Research Use Only

ROD1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01782-1

Capture Antibody Information

Catalog Number:
85061-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242641A6
Reactivity:
human
Immunogen Catalog Number:
Ag5172

Conjugate:
Unconjugated
Full name:
ROD1 regulator of differentiation 1 (S. pombe)
Gene ID:
9991

Detection Antibody Information

Catalog Number:
85061-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242641G5
Reactivity:
human
GenBank:
BC044585
Immunogen Catalog Number:
Ag5172

Conjugate:
Unconjugated
Full name:
ROD1 regulator of differentiation 1 (S. pombe)
Gene ID:
9991

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01782-1 targets ROD1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ROD1 Recombinant antibody, PBS Only (Capture) 85061-1-PBS (242641A6). 100 µg. Concentration 1 mg/ml.

Detection antibody: ROD1 Recombinant antibody, PBS Only (Detector) 85061-3-PBS (242641G5). 100 µg. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01782-1, ROD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85061-1-PBS. Detection antibody: 85061-3-PBS. Standard: Ag5172. Range: 1.563-100 ng/mL.