

For Research Use Only

PML Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01806-2

Capture Antibody Information

Catalog Number:
85048-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242564C12

Reactivity:
human

Immunogen Catalog Number:
Ag15377

Conjugate:
Unconjugated

Full name:
promyelocytic leukemia

Gene ID:
5371

Detection Antibody Information

Catalog Number:
85048-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242564F1

Reactivity:
human

GenBank:
BC020994

Immunogen Catalog Number:
Ag15377

Conjugate:
Unconjugated

Full name:
promyelocytic leukemia

Gene ID:
5371

Applications

Tested Applications:
Cytometric bead array, Sandwich ELISA

Range:
0.781-100 ng/mL (Cytometric Bead Array)
0.156-10 ng/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01806-2 targets PML in immunoassays as a matched antibody pair. Validated in Cytometric bead array, Sandwich ELISA.

Capture antibody: PML Recombinant antibody, PBS Only (Capture) 85048-1-PBS (242564C12). 100 µg. Concentration 1 mg/mL.

Detection antibody: PML Recombinant antibody, PBS Only (Detector) 85048-3-PBS (242564F1). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

