

# DCP1A Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01815-1

## Capture Antibody Information

**Catalog Number:**

85047-2-PBS

**Host:**

Rabbit

**Isotype:**

IgG

**Purification Method:**

Protein A purification

**Clone ID:**

242553B3

**Reactivity:**

human

**Immunogen Catalog Number:**

Ag18027

**Conjugate:**

Unconjugated

**Full name:**

DCP1 decapping enzyme homolog A  
(*S. cerevisiae*)

**Gene ID:**

55802

## Detection Antibody Information

**Catalog Number:**

85047-3-PBS

**Host:**

Rabbit

**Isotype:**

IgG

**Purification Method:**

Protein A purification

**Clone ID:**

242553E12

**Reactivity:**

human

**GenBank:**

BC007439

**Immunogen Catalog Number:**

Ag18027

**Conjugate:**

Unconjugated

**Full name:**

DCP1 decapping enzyme homolog A  
(*S. cerevisiae*)

**Gene ID:**

55802

## Applications

**Tested Applications:**

Cytometric bead array

**Range:**

0.313-10 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01815-1 targets DCP1A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DCP1A Recombinant antibody, PBS Only (Capture) 85047-2-PBS (242553B3). 100 µg. Concentration 1 mg/mL.

Detection antibody: DCP1A Recombinant antibody, PBS Only (Detector) 85047-3-PBS (242553E12). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**

Store at -80°C.

**Storage buffer:**

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MPO1815-1, DCP1A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85047-2-PBS. Detection antibody: 85047-3-PBS. Standard: Ag18027. Range: 0.313-10 ng/mL