

For Research Use Only

TINAGL1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01822-1

Capture Antibody Information

Catalog Number:
85140-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242813D8

Reactivity:
human

Immunogen Catalog Number:
Ag2721

Conjugate:
Unconjugated

Full name:
tubulointerstitial nephritis antigen-like 1

Gene ID:
64129

Detection Antibody Information

Catalog Number:
85140-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242813F10

Reactivity:
human

GenBank:
BC009048

Immunogen Catalog Number:
Ag2721

Conjugate:
Unconjugated

Full name:
tubulointerstitial nephritis antigen-like 1

Gene ID:
64129

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01822-1 targets TINAGL1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TINAGL1 Recombinant antibody, PBS Only (Capture) 85140-1-PBS (242813D8). 100 µg. Concentration 1 mg/ml.

Detection antibody: TINAGL1 Recombinant antibody, PBS Only (Detector) 85140-2-PBS (242813F10). 100 µg. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01822-1, TINAGL1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85140-1-PBS. Detection antibody: 85140-2-PBS. Standard: Ag2721. Range: 1.563-100 ng/mL