

For Research Use Only

# TCL1A Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01835-3

## Capture Antibody Information

Catalog Number:  
85129-4-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
242623E3

Reactivity:  
human

Immunogen Catalog Number:  
Ag0786

Conjugate:  
Unconjugated

Full name:  
T-cell leukemia/lymphoma 1A

Gene ID:  
8115

## Detection Antibody Information

Catalog Number:  
85129-2-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
242623F12

Reactivity:  
human

GenBank:  
BC005831

Immunogen Catalog Number:  
Ag0786

Conjugate:  
Unconjugated

Full name:  
T-cell leukemia/lymphoma 1A

Gene ID:  
8115

## Applications

Tested Applications:  
Sandwich ELISA

Range:  
39.1-2500 pg/mL (Sandwich ELISA)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01835-3 targets TCL1A in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: TCL1A Recombinant antibody, PBS Only (Capture) 85129-4-PBS (242623E3). 100 µg. Concentration 1 mg/mL.

Detection antibody: TCL1A Recombinant antibody, PBS Only (Capture/Detector) 85129-2-PBS (242623F12). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only



For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP01835-3, Human TCL1A Recombinant Matched Antibody Pair - PBS only. 85129-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag0786. 85129-2-PBS was HRP conjugated as the detection antibody. Range: 39.1-2500 pg/mL.

The mean TCL1A concentration was determined to be 627.84 ng/mL in Raji cell extract based on a 3.0 mg/mL extract load, 5,087.48 ng/mL in Ramos cell extract based on a 1.4 mg/mL extract load and 6,616.70 ng/mL in Daudi cell extract based on a 3.6 mg/mL extract load.