

For Research Use Only



FOXA1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01842-3

Capture Antibody Information

Catalog Number:

85125-4-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

242658F4

Reactivity:

human

Immunogen Catalog Number:

Ag14243

Conjugate:

Unconjugated

Full name:

forkhead box A1

Gene ID:

3169

Detection Antibody Information

Catalog Number:

85125-1-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

242658A10

Reactivity:

human

GenBank:

BC033890

Immunogen Catalog Number:

Ag14243

Conjugate:

Unconjugated

Full name:

forkhead box A1

Gene ID:

3169

Applications

Tested Applications:

Sandwich ELISA

Range:

0.195-12.5 ng/mL (Sandwich ELISA)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01842-3 targets FOXA1 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: FOXA1 Recombinant antibody, PBS Only (Capture) 85125-4-PBS (242658F4). 100 µg. Concentration 1 mg/mL.

Detection antibody: FOXA1 Recombinant antibody, PBS Only (Capture/Detector) 85125-1-PBS (242658A10). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:



T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01842-3, Human FOXA1 Recombinant Matched Antibody Pair - PBS only. 85125-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag14243. 85125-1-PBS was HRP conjugated as the detection antibody. Range: 0.195-12.5 ng/mL.

The mean FOXA1 concentration was determined to be 32.58 ng/mL in HepG2 cell extract based on a 1.00 mg/mL extract load, 90.93 ng/mL in MDA-MB-468 cell extract based on a 3.20 mg/mL extract load and 341.89 ng/mL in MCF-7 cell extract based on a 3.50 mg/mL extract load.