

For Research Use Only

Chromogranin A Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01975-2

Capture Antibody Information

Catalog Number:
85605-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
243048H7
Reactivity:
human
Immunogen Catalog Number:
Ag19780

Conjugate:
Unconjugated
Full name:
chromogranin A (parathyroid
secretory protein 1)
Gene ID:
1113

Detection Antibody Information

Catalog Number:
85605-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
243048E8
Reactivity:
human
GenBank:
BC006459
Immunogen Catalog Number:
Ag19780

Conjugate:
Unconjugated
Full name:
chromogranin A (parathyroid
secretory protein 1)
Gene ID:
1113

Applications

Tested Applications:
Cytometric bead array

Range:
3.125-100 ng/mL (Cytometric Bead
Array)

Recommended Dilutions:
It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP01975-2 targets Chromogranin A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CHGA Recombinant antibody, PBS Only (Capture) 85605-2-PBS (243048H7). 100 µg. Concentration 1 mg/mL.

Detection antibody: CHGA Recombinant antibody, PBS Only (Detector) 85605-1-PBS (243048E8). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

**This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.**

Selected Validation Data



Cytometric bead array standard curve of MP01975-2, Chromogranin A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85605-2-PBS. Detection antibody: 85605-1-PBS. Standard: Ag19780. Range: 3.125-100 ng/mL