

Mouse Jam-C Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01981-1**

Capture Antibody Information

Catalog Number:
85615-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242630B1
Reactivity:
mouse

Conjugate:
Unconjugated
Full name:
junction adhesion molecule 3
Gene ID:
83964

Detection Antibody Information

Catalog Number:
85615-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242630E1
Reactivity:
mouse
GenBank:
NM_023277.4

Conjugate:
Unconjugated
Full name:
junction adhesion molecule 3
Gene ID:
83964

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01981-1 targets Jam-C in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse Jam3 Recombinant antibody, PBS Only (Capture) 85615-2-PBS (242630B1). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse Jam3 Recombinant antibody, PBS Only (Detector) 85615-1-PBS (242630E1). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01981-1, MOUSE Jam-C Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85615-2-PBS. Detection antibody: 85615-1-PBS. Standard: Eg2661. Range: 0.625-80 ng/mL.