

For Research Use Only

# ETV5 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50024-1

## Capture Antibody Information

Catalog Number:  
66657-1-PBS

Host:  
Mouse

Isotype:  
IgG2a

Purification Method:  
Protein A purification

Clone ID:  
2A4A4

Reactivity:  
Human

Immunogen Catalog Number:  
Ag3708

Conjugate:  
Unconjugated

Full name:  
ets variant 5

Gene ID:  
2119

## Detection Antibody Information

Catalog Number:  
66657-2-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G purification

Clone ID:  
4E3C2

Reactivity:  
Human

GenBank:  
BC007333

Immunogen Catalog Number:  
Ag3708

Conjugate:  
Unconjugated

Full name:  
ets variant 5

Gene ID:  
2119

## Applications

Tested Applications:  
Sandwich ELISA

Range:  
0.313-20 ng/m (Sandwich ELISA)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50024-1 targets ETV5 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: ETV5 Monoclonal antibody, PBS Only (Capture) 66657-1-PBS (2A4A4). 100 µg. Concentration 1 mg/ml.

Detection antibody: ETV5 Monoclonal antibody, PBS Only (Detector) 66657-2-PBS (4E3C2). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP50024-1, ETV5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66657-1-PBS. Detection antibody: HRP-conjugated 66657-2-PBS. Standard: Ag3708. Range: 0.313-20 ng/m.