

For Research Use Only

TFF1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50027-1

Capture Antibody Information

Catalog Number:
68733-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
2A8D4

Reactivity:
Human

Immunogen Catalog Number:
Ag4541

Conjugate:
Unconjugated

Full name:
trefoil factor 1

Gene ID:
7031

Detection Antibody Information

Catalog Number:
68733-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
2A5D3

Reactivity:
Human

GenBank:
BC032811

Immunogen Catalog Number:
Ag4541

Conjugate:
Unconjugated

Full name:
trefoil factor 1

Gene ID:
7031

Applications

Tested Applications:
Sandwich ELISA

Range:
0.625-20 ng/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50027-1 targets TFF1 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: TFF1 Monoclonal antibody, PBS Only (Capture) 68733-1-PBS (2A8D4). 100 µg. Concentration 1 mg/mL.

Detection antibody: TFF1 Monoclonal antibody, PBS Only (Detector) 68733-2-PBS (2A5D3). 100 µg. Concentration 1 mg/mL.

Alternative TFF1 matched antibody pairs: MP00469-1

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50027-1, TFF1
Monoclonal Matched Antibody Pair, PBS Only.
Capture antibody: 68733-1-PBS. Detection
antibody: HRP-conjugated 68733-2-PBS. Standard:
Ag4541. Range: 0.625-20 ng/mL