

For Research Use Only

TNFSF18 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50040-1

Capture Antibody Information

Catalog Number:
68684-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein A purification

Clone ID:
1G11F11
Reactivity:
Human

Conjugate:
Unconjugated
Full name:
tumor necrosis factor (ligand) superfamily, member 18
Gene ID:
8995

Detection Antibody Information

Catalog Number:
68684-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein A purification

Clone ID:
1C7F5
Reactivity:
Human
GenBank:
BC069319

Conjugate:
Unconjugated
Full name:
tumor necrosis factor (ligand) superfamily, member 18
Gene ID:
8995

Applications

Tested Applications:
Sandwich ELISA

Range:
62.5-4000 pg/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50040-1 targets TNFSF18 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: TNFSF18 Monoclonal antibody, PBS Only (Capture) 68684-1-PBS (1G11F11). 100 µg. Concentration 1 mg/ml.

Detection antibody: TNFSF18 Monoclonal antibody, PBS Only (Detector) 68684-2-PBS (1C7F5). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50040-1, TNFSF18 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68684-1-PBS. Detection antibody: 68684-2-PBS. Standard: Eg32094. Range: 62.5-4000 pg/mL.