

For Research Use Only

TFEB Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50099-2

Capture Antibody Information

Catalog Number:

68632-3-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G purification

Clone ID:

1C7B2

Reactivity:

Human

Immunogen Catalog Number:

Ag24850

Conjugate:

Unconjugated

Full name:

transcription factor EB

Gene ID:

7942

Detection Antibody Information

Catalog Number:

68632-4-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G purification

Clone ID:

1B9C6

Reactivity:

Human

GenBank:

BC032448

Immunogen Catalog Number:

Ag24850

Conjugate:

Unconjugated

Full name:

transcription factor EB

Gene ID:

7942

Applications

Tested Applications:

Cytometric bead array

Range:

0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50099-2 targets TFEB in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TFEB Monoclonal antibody, PBS Only (Capture) 68632-3-PBS (1C7B2). 100 µg. Concentration 1 mg/ml.

Detection antibody: TFEB Monoclonal antibody, PBS Only (Detector) 68632-4-PBS (1B9C6). 100 µg. Concentration 1 mg/ml.

Alternative TFEB matched antibody pairs: MP50099-1

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50099-2, TFEB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68632-3-PBS. Detection antibody: 68632-4-PBS. Standard: Ag24850. Range: 0.781-100 ng/mL.