

# GDF15 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50133-10

## Capture Antibody Information

**Catalog Number:**

68768-6-PBS

**Host:**

Mouse

**Isotype:**

IgG1

**Purification Method:**

Protein G Magarose purification

**Clone ID:**

1F6E2

**Reactivity:**

Human

**Immunogen Catalog Number:**

Ag33750

**Conjugate:**

Unconjugated

**Full name:**

growth differentiation factor 15

**Gene ID:**

9518

## Detection Antibody Information

**Catalog Number:**

68768-2-PBS

**Host:**

Mouse

**Isotype:**

IgG1

**Purification Method:**

Protein G Magarose purification

**Clone ID:**

1B5D12

**Reactivity:**

Human

**GenBank:**

BC008962

**Immunogen Catalog Number:**

Ag33750

**Conjugate:**

Unconjugated

**Full name:**

growth differentiation factor 15

**Gene ID:**

9518

## Applications

**Tested Applications:**

Sandwich ELISA

**Range:**

1.563-100 ng/mL (Sandwich ELISA)

**Recommended Dilutions:**

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50133-10 targets GDF15 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: GDF15 Monoclonal antibody, PBS Only (Capture/Detector) 68768-6-PBS (1F6E2). 100 µg. Concentration 1 mg/ml.

Detection antibody: GDF15 Monoclonal antibody, PBS Only (Capture/Detector) 68768-2-PBS (1B5D12). 100 µg. Concentration 1 mg/ml.

Alternative GDF15 matched antibody pairs: MP50133-1, MP50133-11, MP50133-2, MP50133-3, MP50133-4, MP50133-5, MP50133-6, MP50133-7, MP50133-8, MP50133-9

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**

Store at -80°C.

**Storage buffer:**

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)

W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Sandwich ELISA standard curve of MP50133-10, GDF15 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68768-6-PBS. Detection antibody: HRP-conjugated 68768-2-PBS. Standard: ag33750. Range: 1.563-100 ng/mL.