

For Research Use Only

HMGC51 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50145-4

Capture Antibody Information

Catalog Number:
68779-6-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1E10A8
Reactivity:
Human
Immunogen Catalog Number:
Ag11937

Conjugate:
Unconjugated
Full name:
3-hydroxy-3-methylglutaryl-Coenzyme A synthase 1 (soluble)
Gene ID:
3157

Detection Antibody Information

Catalog Number:
68779-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
2E11C4
Reactivity:
Human
GenBank:
BC083514
Immunogen Catalog Number:
Ag11937

Conjugate:
Unconjugated
Full name:
3-hydroxy-3-methylglutaryl-Coenzyme A synthase 1 (soluble)
Gene ID:
3157

Applications

Tested Applications:
Cytometric bead array, Sandwich ELISA

Range:
0.098-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50145-4 targets HMGC51 in immunoassays as a matched antibody pair. Validated in Cytometric bead array, Sandwich ELISA.

Capture antibody: HMGC51 Monoclonal antibody, PBS Only (Capture) 68779-6-PBS (1E10A8). 100 µg. Concentration 1 mg/mL.

Detection antibody: HMGC51 Monoclonal antibody, PBS Only (Detector) 68779-2-PBS (2E11C4). 100 µg. Concentration 1 mg/mL.

Alternative HMGC51 matched antibody pairs: MP50145-1, MP50145-2, MP50145-3

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50145-4, HMGC51 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68779-6-PBS. Detection antibody: 68779-2-PBS. Standard: Ag11937. Range: 0.098-100 ng/mL.