

For Research Use Only

BSAP,PAX5 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50228-2

Capture Antibody Information

Catalog Number:
60349-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein A purification

Clone ID:
2H2D6

Reactivity:
human

Immunogen Catalog Number:
Ag16023

Conjugate:
Unconjugated

Full name:
paired box 5

Gene ID:
5079

Detection Antibody Information

Catalog Number:
60349-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
5E1E7

Reactivity:
Human

GenBank:
BC156927

Immunogen Catalog Number:
Ag16023

Conjugate:
Unconjugated

Full name:
paired box 5

Gene ID:
5079

Applications

Tested Applications:
Cytometric bead array

Range:
3.125-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50228-2 targets BSAP,PAX5 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BSAP,PAX5 Monoclonal antibody, PBS Only (Capture) 60349-1-PBS (2H2D6). 100 µg. Concentration 1 mg/ml.

Detection antibody: BSAP,PAX5 Monoclonal antibody, PBS Only (Detector) 60349-3-PBS (5E1E7). 100 µg. Concentration 1 mg/ml.

Alternative BSAP,PAX5 matched antibody pairs: MP50228-1, MP50228-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50228-2, BSAP,PAX5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60349-1-PBS. Detection antibody: 60349-3-PBS. Standard:Ag16023. Range: 3.125-100 ng/mL