

For Research Use Only

IL4R Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50242-2

Capture Antibody Information

Catalog Number:

67051-4-PBS

Host:

Mouse

Isotype:

IgM

Purification Method:

Euglobulin precipitation

Clone ID:

1B10E9

Reactivity:

Human

Immunogen Catalog Number:

Ag28411

Conjugate:

Unconjugated

Full name:

interleukin 4 receptor

Gene ID:

3566

Detection Antibody Information

Catalog Number:

67051-3-PBS

Host:

Mouse

Isotype:

IgM

Purification Method:

Euglobulin precipitation

Clone ID:

3F12E8

Reactivity:

Human

GenBank:

NM_000418

Immunogen Catalog Number:

Ag28411

Conjugate:

Unconjugated

Full name:

interleukin 4 receptor

Gene ID:

3566

Applications

Tested Applications:

N/A

Range:

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50242-2 targets IL4R in immunoassays as a matched antibody pair. Validated in N/A.

Capture antibody: IL4R Monoclonal antibody, PBS Only (Capture) 67051-4-PBS (1B10E9). 100 µg. Concentration 1 mg/ml.

Detection antibody: IL4R Monoclonal antibody, PBS Only (Capture/Detector) 67051-3-PBS (3F12E8). 100 µg. Concentration 1 mg/ml.

Alternative IL4R matched antibody pairs: MP50242-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50242-2, IL4R Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67051-4-PBS. Detection antibody: 67051-3-PBS. Standard: Ag28411. Range: 0.781-100 ng/mL.