

For Research Use Only

CHMP1A Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50256-1

Capture Antibody Information

Catalog Number:
68862-1-PBS
Host:
Mouse
Isotype:
IgG3
Purification Method:
Protein A Magarose purification

Clone ID:
1H9D10
Reactivity:
Human
Immunogen Catalog Number:
Ag8376

Conjugate:
Unconjugated
Full name:
chromatin modifying protein 1A
Gene ID:
5119

Detection Antibody Information

Catalog Number:
68862-2-PBS
Host:
Mouse
Isotype:
IgG2a
Purification Method:
Protein A Magarose purification

Clone ID:
1D7E7
Reactivity:
Human
GenBank:
BC007527
Immunogen Catalog Number:
Ag8376

Conjugate:
Unconjugated
Full name:
chromatin modifying protein 1A
Gene ID:
5119

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50256-1 targets CHMP1A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CHMP1A Monoclonal antibody, PBS Only (Capture) 68862-1-PBS (1H9D10). 100 µg. Concentration 1 mg/ml.

Detection antibody: CHMP1A Monoclonal antibody, PBS Only (Detector) 68862-2-PBS (1D7E7). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50256-1, CHMP1A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68862-1-PBS. Detection antibody: 68862-2-PBS. Standard: Ag8376. Range: 0.781-100 ng/mL.