

For Research Use Only

p504S,AMACR Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50263-1

Capture Antibody Information

Catalog Number:
60240-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
5F4D7
Reactivity:
human
Immunogen Catalog Number:
Ag8720

Conjugate:
Unconjugated
Full name:
alpha-methylacyl-CoA racemase
Gene ID:
23600

Detection Antibody Information

Catalog Number:
60240-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
5F5D6
Reactivity:
human
GenBank:
BC009471
Immunogen Catalog Number:
Ag8720

Conjugate:
Unconjugated
Full name:
alpha-methylacyl-CoA racemase
Gene ID:
23600

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50263-1 targets p504S,AMACR in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: p504S,AMACR Monoclonal antibody, PBS Only (Capture) 60240-2-PBS (5F4D7). 100 µg. Concentration 1 mg/ml.
Detection antibody: p504S,AMACR Monoclonal antibody, PBS Only (Capture/Detector) 60240-3-PBS (5F5D6). 100 µg. Concentration 1 mg/ml.
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50263-1, p504S,AMACR Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60240-2-PBS. Detection antibody: 60240-3-PBS. Standard:Ag8720. Range: 1.563-200 ng/mL