

For Research Use Only

SPHK2 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50264-3

Capture Antibody Information

Catalog Number:
68872-2-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A purification

Clone ID:
3F7H5

Reactivity:
human

Immunogen Catalog Number:
Ag10569

Conjugate:
Unconjugated

Full name:
sphingosine kinase 2

Gene ID:
56848

Detection Antibody Information

Catalog Number:
68871-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
1B2H4

Reactivity:
human

GenBank:
BC010671

Immunogen Catalog Number:
Ag33745

Conjugate:
Unconjugated

Full name:
sphingosine kinase 2

Gene ID:
56848

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-25 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50264-3 targets SPHK2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SPHK2 Monoclonal antibody, PBS Only (Capture) 68872-2-PBS (3F7H5). 100 µg. Concentration 1 mg/mL.

Detection antibody: SPHK2 Monoclonal antibody, PBS Only (Detector) 68871-3-PBS (1B2H4). 100 µg. Concentration 1 mg/mL.

Alternative SPHK2 matched antibody pairs: MP50264-1, MP50264-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50264-3, SPHK2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68872-2-PBS. Detection antibody: 68871-3-PBS. Standard: Ag10569. Range: 0.098-25 ng/mL