

For Research Use Only

STK32A Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50267-2

Capture Antibody Information

Catalog Number:
67528-1-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A purification

Clone ID:
1H2A6

Reactivity:
human

Immunogen Catalog Number:
Ag10942

Conjugate:
Unconjugated

Full name:
serine/threonine kinase 32A

Gene ID:
202374

Detection Antibody Information

Catalog Number:
67528-3-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A Magarose purification

Clone ID:
2B6E5

Reactivity:
human

GenBank:
BC021666

Immunogen Catalog Number:
Ag10942

Conjugate:
Unconjugated

Full name:
serine/threonine kinase 32A

Gene ID:
202374

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50267-2 targets STK32A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: STK32A Monoclonal antibody, PBS Only (Capture) 67528-1-PBS (1H2A6). 100 µg. Concentration 1 mg/ml.

Detection antibody: STK32A Monoclonal antibody, PBS Only (Detector) 67528-3-PBS (2B6E5). 100 µg. Concentration 1 mg/ml.

Alternative STK32A matched antibody pairs: MP50267-1, MP50267-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50267-2, STK32A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67528-1-PBS. Detection antibody: 67528-3-PBS. Standard: Ag10942. Range: 1.563-200 ng/mL.