

DLL4 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50281-2

Capture Antibody Information

Catalog Number:

68882-1-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G Magarose purification

Clone ID:

2A3H2

Reactivity:

human

Immunogen Catalog Number:

Ag17483

Conjugate:

Unconjugated

Full name:

delta-like 4 (Drosophila)

Gene ID:

54567

Detection Antibody Information

Catalog Number:

68882-3-PBS

Host:

Mouse

Isotype:

IgG2a

Purification Method:

Protein A Magarose purification

Clone ID:

1F3D12

Reactivity:

human

GenBank:

BC106950

Immunogen Catalog Number:

Ag17483

Conjugate:

Unconjugated

Full name:

delta-like 4 (Drosophila)

Gene ID:

54567

Applications

Tested Applications:

N/A

Range:

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50281-2 targets DLL4 in immunoassays as a matched antibody pair. Validated in N/A.

Capture antibody: DLL4 Monoclonal antibody, PBS Only (Capture) 68882-1-PBS (2A3H2). 100 µg. Concentration 1 mg/ml.

Detection antibody: DLL4 Monoclonal antibody, PBS Only (Detector) 68882-3-PBS (1F3D12). 100 µg. Concentration 1 mg/ml.

Alternative DLL4 matched antibody pairs: MP50281-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50281-2, DLL4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68882-1-PBS. Detection antibody: 68882-3-PBS. Standard: Ag16140. Range: 0.391-12.5 ng/mL.