

For Research Use Only

TRIM11 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50306-1

Capture Antibody Information

Catalog Number:
68904-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
2D6D1
Reactivity:
human
Immunogen Catalog Number:
Ag1308

Conjugate:
Unconjugated
Full name:
tripartite motif-containing 11
Gene ID:
81559

Detection Antibody Information

Catalog Number:
68904-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
2G5E10
Reactivity:
human
GenBank:
BC011629
Immunogen Catalog Number:
Ag1308

Conjugate:
Unconjugated
Full name:
tripartite motif-containing 11
Gene ID:
81559

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50306-1 targets TRIM11 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TRIM11 Monoclonal antibody, PBS Only (Capture) 68904-1-PBS (2D6D1). 100 µg. Concentration 1 mg/ml.

Detection antibody: TRIM11 Monoclonal antibody, PBS Only (Detector) 68904-2-PBS (2G5E10). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50306-1, TRIM11 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68904-1-PBS. Detection antibody: 68904-2-PBS. Standard:Ag1308. Range: 0.781-50 ng/mL