

NPC2 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50350-2

Capture Antibody Information

Catalog Number:

68927-3-PBS

Host:

Mouse

Isotype:

IgG2b

Purification Method:

Protein A Magarose purification

Clone ID:

2C7E3

Reactivity:

human

Immunogen Catalog Number:

Ag34710

Conjugate:

Unconjugated

Full name:

Niemann-Pick disease, type C2

Gene ID:

10577

Detection Antibody Information

Catalog Number:

68927-4-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G Magarose purification

Clone ID:

2D8D11

Reactivity:

human

GenBank:

BC002532

Immunogen Catalog Number:

Ag34710

Conjugate:

Unconjugated

Full name:

Niemann-Pick disease, type C2

Gene ID:

10577

Applications

Tested Applications:

Cytometric bead array

Range:

1.563-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50350-2 targets NPC2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NPC2 Monoclonal antibody, PBS Only (Capture) 68927-3-PBS (2C7E3). 100 µg. Concentration 1 mg/ml.

Detection antibody: NPC2 Monoclonal antibody, PBS Only (Detector) 68927-4-PBS (2D8D11). 100 µg. Concentration 1 mg/ml.

Alternative NPC2 matched antibody pairs: MP50350-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50350-2, NPC2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68927-3-PBS. Detection antibody: 68927-4-PBS. Standard: Ag34710. Range: 1.563-100 ng/mL.