

For Research Use Only

TRIM21 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50371-4

Capture Antibody Information

Catalog Number:

67136-1-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein A purification

Clone ID:

1E8B9

Reactivity:

human, mouse, rat

Immunogen Catalog Number:

Ag28377

Conjugate:

Unconjugated

Full name:

tripartite motif-containing 21

Gene ID:

6737

Detection Antibody Information

Catalog Number:

67136-5-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G purification

Clone ID:

1F4H7

Reactivity:

human

GenBank:

BC010861

Immunogen Catalog Number:

Ag28377

Conjugate:

Unconjugated

Full name:

tripartite motif-containing 21

Gene ID:

6737

Applications

Tested Applications:

Cytometric bead array

Range:

0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50371-4 targets TRIM21 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TRIM21 Monoclonal antibody, PBS Only (Capture) 67136-1-PBS (1E8B9). 100 µg. Concentration 1 mg/ml.

Detection antibody: TRIM21 Monoclonal antibody, PBS Only (Detector) 67136-5-PBS (1F4H7). 100 µg. Concentration 1 mg/ml.

Alternative TRIM21 matched antibody pairs: MP00505-1, MP50371-1, MP50371-2, MP50371-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50371-4, TRIM21 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67136-1-PBS. Detection antibody: 67136-5-PBS. Standard: Ag28377. Range: 0.781-100 ng/mL.