

For Research Use Only

PRKCQ Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50373-2

Capture Antibody Information

Catalog Number:
68939-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
1E2E2

Reactivity:
human

Immunogen Catalog Number:
Ag26331

Conjugate:
Unconjugated

Full name:
protein kinase C, theta

Gene ID:
5588

Detection Antibody Information

Catalog Number:
68939-4-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
3E2C3

Reactivity:
human

GenBank:
BC101465

Immunogen Catalog Number:
Ag26331

Conjugate:
Unconjugated

Full name:
protein kinase C, theta

Gene ID:
5588

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-25 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50373-2 targets PRKCQ in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PRKCQ Monoclonal antibody, PBS Only (Capture) 68939-3-PBS (1E2E2). 100 µg. Concentration 1 mg/mL.

Detection antibody: PRKCQ Monoclonal antibody, PBS Only (Detector) 68939-4-PBS (3E2C3). 100 µg. Concentration 1 mg/mL.

Alternative PRKCQ matched antibody pairs: MP50373-1, MP50373-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50373-2, PRKCC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68939-3-PBS. Detection antibody: 68939-4-PBS. Standard: Ag26331. Range: 0.781-25 ng/mL