

For Research Use Only

BACH1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50384-1

Capture Antibody Information

Catalog Number:
66762-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein A purification

Clone ID:
2D12A4

Reactivity:
human, mouse, rat, pig, rabbit

Immunogen Catalog Number:
Ag6011

Conjugate:
Unconjugated

Full name:
BTB and CNC homology 1, basic
leucine zipper transcription factor 1

Gene ID:
571

Detection Antibody Information

Catalog Number:
66762-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
4E4A5

Reactivity:
human

GenBank:
BC063307

Immunogen Catalog Number:
Ag6011

Conjugate:
Unconjugated

Full name:
BTB and CNC homology 1, basic
leucine zipper transcription factor 1

Gene ID:
571

Applications

Tested Applications:
Cytometric bead array

Range:
6.25-200 ng/mL (Cytometric Bead
Array)

Recommended Dilutions:
It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP50384-1 targets BACH1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BACH1 Monoclonal antibody, PBS Only (Capture) 66762-1-PBS (2D12A4). 100 µg. Concentration 1 mg/ml.

Detection antibody: BACH1 Monoclonal antibody, PBS Only (Capture/Detector) 66762-2-PBS (4E4A5). 100 µg. Concentration 1 mg/ml.

Alternative BACH1 matched antibody pairs: MP50384-2, MP80007-1, MP80007-2, MP80007-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50384-1, BACH1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66762-1-PBS. Detection antibody: 66762-2-PBS. Standard: Ag6011. Range: 6.25-200 ng/mL