

For Research Use Only

# ADAM17 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50392-3

## Capture Antibody Information

Catalog Number:

68725-2-PBS

Host:

Mouse

Isotype:

IgG2a

Purification Method:

Protein A purification

Clone ID:

1D5A9

Reactivity:

human

Immunogen Catalog Number:

Ag34022

Conjugate:

Unconjugated

Full name:

ADAM metalloproteinase domain 17

Gene ID:

6868

## Detection Antibody Information

Catalog Number:

68725-4-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G purification

Clone ID:

1E1G7

Reactivity:

human

GenBank:

BC136783

Immunogen Catalog Number:

Ag34022

Conjugate:

Unconjugated

Full name:

ADAM metalloproteinase domain 17

Gene ID:

6868

## Applications

Tested Applications:

Cytometric bead array

Range:

0.098-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50392-3 targets ADAM17 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ADAM17 Monoclonal antibody, PBS Only (Capture) 68725-2-PBS (1D5A9). 100 µg. Concentration 1 mg/ml.

Detection antibody: ADAM17 Monoclonal antibody, PBS Only (Detector) 68725-4-PBS (1E1G7). 100 µg. Concentration 1 mg/ml.

Alternative ADAM17 matched antibody pairs: MP50392-1, MP50392-2, MP50392-4

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)

W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50392-3, ADAM17 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68725-2-PBS. Detection antibody: 68725-4-PBS. Standard: Ag34022. Range: 0.098-100 ng/mL.