

For Research Use Only

# PAX8 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50408-6

## Capture Antibody Information

Catalog Number:  
60144-3-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G Magarose purification

Clone ID:  
3C2G10

Reactivity:  
human

Immunogen Catalog Number:  
Ag0306

Conjugate:  
Unconjugated

Full name:  
paired box 8

Gene ID:  
7849

## Detection Antibody Information

Catalog Number:  
68979-1-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G Magarose purification

Clone ID:  
1F4C3

Reactivity:  
human

GenBank:  
BC001060

Immunogen Catalog Number:  
Ag16019

Conjugate:  
Unconjugated

Full name:  
paired box 8

Gene ID:  
7849

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.781-25 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50408-6 targets PAX8 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PAX8 Monoclonal antibody, PBS Only (Capture) 60144-3-PBS (3C2G10). 100 µg. Concentration 1 mg/ml.

Detection antibody: PAX8 Monoclonal antibody, PBS Only (Detector) 68979-1-PBS (1F4C3). 100 µg. Concentration 1 mg/ml.

Alternative PAX8 matched antibody pairs: MP50408-1, MP50408-2, MP50408-3, MP50408-4, MP50408-5

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50408-6, PAX8 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60144-3-PBS. Detection antibody: 68979-1-PBS. Standard: Ag0306. Range: 0.781-25 ng/mL