

For Research Use Only

SLC25A37 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50411-1

Capture Antibody Information

Catalog Number:
67593-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
2A3E3

Reactivity:
human

Immunogen Catalog Number:
Ag26216

Conjugate:
Unconjugated

Full name:
solute carrier family 25, member 37

Gene ID:
51312

Detection Antibody Information

Catalog Number:
67593-3-PBS

Host:
Mouse

Isotype:
IgG3

Purification Method:
Protein A Magarose purification

Clone ID:
1H2B2

Reactivity:
human

GenBank:
BC132799

Immunogen Catalog Number:
Ag26216

Conjugate:
Unconjugated

Full name:
solute carrier family 25, member 37

Gene ID:
51312

Applications

Tested Applications:
Cytometric bead array

Range:
6.25-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50411-1 targets SLC25A37 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SLC25A37 Monoclonal antibody, PBS Only (Capture) 67593-2-PBS (2A3E3). 100 µg. Concentration 1 mg/ml.

Detection antibody: SLC25A37 Monoclonal antibody, PBS Only (Detector) 67593-3-PBS (1H2B2). 100 µg. Concentration 1 mg/ml.

Alternative SLC25A37 matched antibody pairs: MP50411-2, MP50411-3, MP50411-4

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50411-1, SLC25A37 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67593-2-PBS. Detection antibody: 67593-3-PBS. Standard: Ag26216. Range: 6.25-200 ng/mL