

For Research Use Only

# MYH1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50449-2

## Capture Antibody Information

**Catalog Number:**  
67299-4-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G Magarose purification

**Clone ID:**  
3B3A5  
**Reactivity:**  
human  
**Immunogen Catalog Number:**  
Ag17129

**Conjugate:**  
Unconjugated  
**Full name:**  
myosin, heavy chain 1, skeletal muscle, adult  
**Gene ID:**  
4619

## Detection Antibody Information

**Catalog Number:**  
67299-5-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G purification

**Clone ID:**  
1E5C4  
**Reactivity:**  
human  
**GenBank:**  
BC114545  
**Immunogen Catalog Number:**  
Ag17129

**Conjugate:**  
Unconjugated  
**Full name:**  
myosin, heavy chain 1, skeletal muscle, adult  
**Gene ID:**  
4619

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
3.125-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50449-2 targets MYH1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.  
Capture antibody: MYH1 Monoclonal antibody, PBS Only (Capture) 67299-4-PBS (3B3A5). 100 µg. Concentration 1 mg/ml.  
Detection antibody: MYH1 Monoclonal antibody, PBS Only (Detector) 67299-5-PBS (1E5C4). 100 µg. Concentration 1 mg/ml.  
Alternative MYH1 matched antibody pairs: MP50449-1  
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.  
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.  
Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50449-2, MYH1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67299-4-PBS. Detection antibody: 67299-5-PBS. Standard: Ag17129. Range: 3.125-100 ng/mL.