

For Research Use Only

# NUP210 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50488-1

## Capture Antibody Information

<b>Catalog Number:</b> 68996-1-PBS	<b>Clone ID:</b> 1C6A1	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> nucleoporin 210kDa
<b>Isotype:</b> IgG2a	<b>Immunogen Catalog Number:</b> Ag10494	<b>Gene ID:</b> 23225
<b>Purification Method:</b> Protein A Magarose purification		

## Detection Antibody Information

<b>Catalog Number:</b> 68996-2-PBS	<b>Clone ID:</b> 2D8D5	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> nucleoporin 210kDa
<b>Isotype:</b> IgG1	<b>GenBank:</b> BC067089	<b>Gene ID:</b> 23225
<b>Purification Method:</b> Protein G Magarose purification		
<b>Immunogen Catalog Number:</b> Ag10494		

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 0.391-100 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

## Product Information

MP50488-1 targets NUP210 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NUP210 Monoclonal antibody, PBS Only (Capture) 68996-1-PBS (1C6A1). 100 µg. Concentration 1 mg/ml.

Detection antibody: NUP210 Monoclonal antibody, PBS Only (Detector) 68996-2-PBS (2D8D5). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50488-1, NUP210 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68996-1-PBS. Detection antibody: 68996-2-PBS. Standard: Ag10494. Range: 0.391-100 ng/mL.