

For Research Use Only

CHCHD4 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50491-2

Capture Antibody Information

Catalog Number:
66718-4-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1G2A11
Reactivity:
human
Immunogen Catalog Number:
Ag15339

Conjugate:
Unconjugated
Full name:
coiled-coil-helix-coiled-coil-helix domain containing 4
Gene ID:
131474

Detection Antibody Information

Catalog Number:
66718-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
3D8B9
Reactivity:
human
GenBank:
BC033775
Immunogen Catalog Number:
Ag15339

Conjugate:
Unconjugated
Full name:
coiled-coil-helix-coiled-coil-helix domain containing 4
Gene ID:
131474

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50491-2 targets CHCHD4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: CHCHD4 Monoclonal antibody, PBS Only (Capture) 66718-4-PBS (1G2A11). 100 µg. Concentration 1 mg/ml.
Detection antibody: CHCHD4 Monoclonal antibody, PBS Only (Detector) 66718-3-PBS (3D8B9). 100 µg. Concentration 1 mg/ml.
Alternative CHCHD4 matched antibody pairs: MP50491-1, MP50491-3
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50491-2, CHCHD4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66718-4-PBS. Detection antibody: 66718-3-PBS. Standard: Ag15339. Range: 0.781-50 ng/mL