

For Research Use Only

# SLC10A3 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50537-1

## Capture Antibody Information

Catalog Number:

60419-1-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G Magarose purification

Clone ID:

1H1A11

Reactivity:

human

Immunogen Catalog Number:

Ag14038

Conjugate:

Unconjugated

Full name:

solute carrier family 10 (sodium/bile acid cotransporter family), member 3

Gene ID:

8273

## Detection Antibody Information

Catalog Number:

60419-2-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G Magarose purification

Clone ID:

4A10B4

Reactivity:

human

GenBank:

BC004966

Immunogen Catalog Number:

Ag14038

Conjugate:

Unconjugated

Full name:

solute carrier family 10 (sodium/bile acid cotransporter family), member 3

Gene ID:

8273

## Applications

Tested Applications:

Cytometric bead array

Range:

0.781-12.5 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50537-1 targets SLC10A3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SLC10A3 Monoclonal antibody, PBS Only (Capture) 60419-1-PBS (1H1A11). 100 µg. Concentration 1 mg/ml.

Detection antibody: SLC10A3 Monoclonal antibody, PBS Only (Detector) 60419-2-PBS (4A10B4). 100 µg. Concentration 1 mg/ml.

Alternative SLC10A3 matched antibody pairs: MP50537-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)

W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50537-1, SLC10A3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60419-1-PBS. Detection antibody: 60419-2-PBS. Standard: Ag14038. Range: 0.781-12.5 ng/mL.