

For Research Use Only

CCDC56 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50539-1

Capture Antibody Information

Catalog Number:
60420-1-PBS

Host:
Mouse

Isotype:
IgG2b

Purification Method:
Protein A purification

Clone ID:
1E10H10

Reactivity:
human

Immunogen Catalog Number:
Ag14730

Conjugate:
Unconjugated

Full name:
coiled-coil domain containing 56

Gene ID:
28958

Detection Antibody Information

Catalog Number:
60420-2-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A purification

Clone ID:
1D8H5

Reactivity:
human

GenBank:
BC002698

Immunogen Catalog Number:
Ag14730

Conjugate:
Unconjugated

Full name:
coiled-coil domain containing 56

Gene ID:
28958

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50539-1 targets CCDC56 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CCDC56 Monoclonal antibody, PBS Only (Capture) 60420-1-PBS (1E10H10). 100 µg. Concentration 1 mg/mL.

Detection antibody: CCDC56 Monoclonal antibody, PBS Only (Detector) 60420-2-PBS (1D8H5). 100 µg. Concentration 1 mg/mL.

Alternative CCDC56 matched antibody pairs: MP50539-2, MP50539-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50539-1, CCDC56 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60420-1-PBS. Detection antibody: 60420-2-PBS. Standard:Ag14730. Range: 0.781-100 ng/mL.