

For Research Use Only

SDHB Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50553-1

Capture Antibody Information

Catalog Number:
67600-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
2C5A5
Reactivity:
human
Immunogen Catalog Number:
Ag29868

Conjugate:
Unconjugated
Full name:
succinate dehydrogenase complex, subunit B, iron sulfur (Ip)
Gene ID:
6390

Detection Antibody Information

Catalog Number:
67600-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
2E7E12
Reactivity:
human
GenBank:
BC007840
Immunogen Catalog Number:
Ag29868

Conjugate:
Unconjugated
Full name:
succinate dehydrogenase complex, subunit B, iron sulfur (Ip)
Gene ID:
6390

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50553-1 targets SDHB in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: SDHB Monoclonal antibody, PBS Only (Capture/Detector) 67600-2-PBS (2C5A5). 100 µg. Concentration 1 mg/ml.
Detection antibody: SDHB Monoclonal antibody, PBS Only (Detector) 67600-3-PBS (2E7E12). 100 µg. Concentration 1 mg/ml.
Alternative SDHB matched antibody pairs: MP50553-2, MP50553-3
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50553-1, SDHB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67600-2-PBS. Detection antibody: 67600-3-PBS. Standard: Ag29868. Range: 1.563-100 ng/mL.