

For Research Use Only

Tim23 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50554-3

Capture Antibody Information

Catalog Number:
67535-5-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1A10E7
Reactivity:
human
Immunogen Catalog Number:
Ag27069

Conjugate:
Unconjugated
Full name:
translocase of inner mitochondrial membrane 23 homolog (yeast)
Gene ID:
10431

Detection Antibody Information

Catalog Number:
67535-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
2E7G5
Reactivity:
human
GenBank:
BC062707
Immunogen Catalog Number:
Ag27069

Conjugate:
Unconjugated
Full name:
translocase of inner mitochondrial membrane 23 homolog (yeast)
Gene ID:
10431

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50554-3 targets Tim23 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: Tim23 Monoclonal antibody, PBS Only (Capture) 67535-5-PBS (1A10E7). 100 µg. Concentration 1 mg/ml.
Detection antibody: Tim23 Monoclonal antibody, PBS Only (Capture/Detector) 67535-3-PBS (2E7G5). 100 µg. Concentration 1 mg/ml.
Alternative Tim23 matched antibody pairs: MP50554-1, MP50554-2
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50554-3, Tim23 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67535-5-PBS. Detection antibody: 67535-3-PBS. Standard: Ag27069. Range: 1.563-100 ng/mL.