

For Research Use Only

VAMP8 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50557-2

Capture Antibody Information

Catalog Number:
60425-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
3B6G5

Reactivity:
human

Immunogen Catalog Number:
Ag8381

Conjugate:
Unconjugated

Full name:
vesicle-associated membrane protein 8 (endobrevin)

Gene ID:
8673

Detection Antibody Information

Catalog Number:
60425-4-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1H9E5

Reactivity:
human

GenBank:
BC001634

Immunogen Catalog Number:
Ag8381

Conjugate:
Unconjugated

Full name:
vesicle-associated membrane protein 8 (endobrevin)

Gene ID:
8673

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50557-2 targets VAMP8 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: VAMP8 Monoclonal antibody, PBS Only (Capture) 60425-3-PBS (3B6G5). 100 µg. Concentration 1 mg/ml.

Detection antibody: VAMP8 Monoclonal antibody, PBS Only (Detector) 60425-4-PBS (1H9E5). 100 µg. Concentration 1 mg/ml.

Alternative VAMP8 matched antibody pairs: MP50557-1, MP50557-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50557-2, VAMP8 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60425-3-PBS. Detection antibody: 60425-4-PBS. Standard: Ag8381. Range: 1.563-100 ng/mL.