

For Research Use Only

# CD63 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50558-2

## Capture Antibody Information

Catalog Number:  
67605-4-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G purification

Clone ID:  
1E12F6

Reactivity:  
human

Immunogen Catalog Number:  
Ag19690

Conjugate:  
Unconjugated

Full name:  
CD63 molecule

Gene ID:  
967

## Detection Antibody Information

Catalog Number:  
67605-3-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G Magarose purification

Clone ID:  
1B2B2

Reactivity:  
human

GenBank:  
BC002349

Immunogen Catalog Number:  
Ag19690

Conjugate:  
Unconjugated

Full name:  
CD63 molecule

Gene ID:  
967

## Applications

Tested Applications:  
Cytometric bead array

Range:  
1.563-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50558-2 targets CD63 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CD63 Monoclonal antibody, PBS Only (Capture) 67605-4-PBS (1E12F6). 100 µg. Concentration 1 mg/ml.

Detection antibody: CD63 Monoclonal antibody, PBS Only (Detector) 67605-3-PBS (1B2B2). 100 µg. Concentration 1 mg/ml.

Alternative CD63 matched antibody pairs: MP50558-1, MP50558-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50558-2, CD63 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67605-4-PBS. Detection antibody: 67605-3-PBS. Standard:Ag19690. Range: 1.563-100 ng/mL.