

For Research Use Only

CPVL Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50562-3

Capture Antibody Information

Catalog Number:
60429-4-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A Magarose purification

Clone ID:
1F3H4

Reactivity:
human

Immunogen Catalog Number:
Ag36547

Conjugate:
Unconjugated

Full name:
carboxypeptidase, vitellogenic-like

Gene ID:
54504

Detection Antibody Information

Catalog Number:
60429-5-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2E10G6

Reactivity:
human

GenBank:
BC016838

Immunogen Catalog Number:
Ag36547

Conjugate:
Unconjugated

Full name:
carboxypeptidase, vitellogenic-like

Gene ID:
54504

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50562-3 targets CPVL in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CPVL Monoclonal antibody, PBS Only (Capture) 60429-4-PBS (1F3H4). 100 µg. Concentration 1 mg/ml.

Detection antibody: CPVL Monoclonal antibody, PBS Only (Detector) 60429-5-PBS (2E10G6). 100 µg. Concentration 1 mg/ml.

Alternative CPVL matched antibody pairs: MP50562-1, MP50562-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50562-3, CPVL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60429-4-PBS. Detection antibody: 60429-5-PBS. Standard: Ag36547. Range: 0.781-100 ng/mL.