

For Research Use Only

ETF A Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50564-1

Capture Antibody Information

Catalog Number:
67661-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1G6E2
Reactivity:
human
Immunogen Catalog Number:
Ag30257

Conjugate:
Unconjugated
Full name:
electron-transfer-flavoprotein, alpha polypeptide
Gene ID:
2108

Detection Antibody Information

Catalog Number:
67661-3-PBS
Host:
Mouse
Isotype:
IgG2b
Purification Method:
Protein A Magarose purification

Clone ID:
1F10C3
Reactivity:
human
GenBank:
BC015526
Immunogen Catalog Number:
Ag30257

Conjugate:
Unconjugated
Full name:
electron-transfer-flavoprotein, alpha polypeptide
Gene ID:
2108

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50564-1 targets ETF A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: ETF A Monoclonal antibody, PBS Only (Capture) 67661-2-PBS (1G6E2). 100 µg. Concentration 1 mg/ml.
Detection antibody: ETF A Monoclonal antibody, PBS Only (Detector) 67661-3-PBS (1F10C3). 100 µg. Concentration 1 mg/ml.
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50564-1, ETFA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67661-2-PBS. Detection antibody: 67661-3-PBS. Standard: Ag30257. Range: 1.563-100 ng/mL.