

For Research Use Only

# OMA1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50571-1

## Capture Antibody Information

**Catalog Number:**  
66572-1-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G purification

**Clone ID:**  
2E11D7  
**Reactivity:**  
human  
**Immunogen Catalog Number:**  
Ag11201

**Conjugate:**  
Unconjugated  
**Full name:**  
OMA1 homolog, zinc metalloproteinase (*S. cerevisiae*)  
**Gene ID:**  
115209

## Detection Antibody Information

**Catalog Number:**  
66572-2-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G purification

**Clone ID:**  
1G12H1  
**Reactivity:**  
human  
**GenBank:**  
BC012915  
**Immunogen Catalog Number:**  
Ag11201

**Conjugate:**  
Unconjugated  
**Full name:**  
OMA1 homolog, zinc metalloproteinase (*S. cerevisiae*)  
**Gene ID:**  
115209

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.781-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50571-1 targets OMA1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: OMA1 Monoclonal antibody, PBS Only (Capture) 66572-1-PBS (2E11D7). 100 µg. Concentration 1 mg/ml.

Detection antibody: OMA1 Monoclonal antibody, PBS Only (Detector) 66572-2-PBS (1G12H1). 100 µg. Concentration 1 mg/ml.

Alternative OMA1 matched antibody pairs: MP50571-2, MP50571-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50571-1, OMA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66572-1-PBS. Detection antibody: 66572-2-PBS. Standard: Ag11201. Range: 0.781-100 ng/mL.