

For Research Use Only

ALCAM Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50574-4

Capture Antibody Information

Catalog Number:
67768-5-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
3C1D12
Reactivity:
human
Immunogen Catalog Number:
Ag16892

Conjugate:
Unconjugated
Full name:
activated leukocyte cell adhesion molecule
Gene ID:
214

Detection Antibody Information

Catalog Number:
67768-7-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2F6D6
Reactivity:
human
GenBank:
BC137097
Immunogen Catalog Number:
Ag16892

Conjugate:
Unconjugated
Full name:
activated leukocyte cell adhesion molecule
Gene ID:
214

Applications

Tested Applications:
Cytometric bead array

Range:
25-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50574-4 targets ALCAM in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: ALCAM Monoclonal antibody, PBS Only (Capture/Detector) 67768-5-PBS (3C1D12). 100 µg. Concentration 1 mg/ml.
Detection antibody: ALCAM Monoclonal antibody, PBS Only (Detector) 67768-7-PBS (2F6D6). 100 µg. Concentration 1 mg/ml.
Alternative ALCAM matched antibody pairs: MP50574-1, MP50574-2, MP50574-3
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50574-4, ALCAM Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67768-5-PBS. Detection antibody: 67768-7-PBS. Standard: Ag16892. Range: 25-200 ng/mL