

ATP5L Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50581-2

Capture Antibody Information

Catalog Number:

60434-1-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G Magarose purification

Clone ID:

2D10F9

Reactivity:

human

Immunogen Catalog Number:

Ag9287

Conjugate:

Unconjugated

Full name:

ATP synthase, H⁺ transporting, mitochondrial FO complex, subunit G

Gene ID:

10632

Detection Antibody Information

Catalog Number:

60434-3-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G Magarose purification

Clone ID:

2H8H3

Reactivity:

human

GenBank:

BC015128

Immunogen Catalog Number:

Ag9287

Conjugate:

Unconjugated

Full name:

ATP synthase, H⁺ transporting, mitochondrial FO complex, subunit G

Gene ID:

10632

Applications

Tested Applications:

Cytometric bead array

Range:

0.391-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50581-2 targets ATP5L in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ATP5L Monoclonal antibody, PBS Only (Capture) 60434-1-PBS (2D10F9). 100 µg. Concentration 1 mg/ml.

Detection antibody: ATP5L Monoclonal antibody, PBS Only (Detector) 60434-3-PBS (2H8H3). 100 µg. Concentration 1 mg/ml.

Alternative ATP5L matched antibody pairs: MP50581-1, MP50581-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50581-2, ATP5L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60434-1-PBS. Detection antibody: 60434-3-PBS. Standard: Ag9287. Range: 0.391-100 ng/mL.