

For Research Use Only

CYP11B1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50582-3

Capture Antibody Information

Catalog Number:
60435-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1C8A4
Reactivity:
human
Immunogen Catalog Number:
Ag17158

Conjugate:
Unconjugated
Full name:
cytochrome P450, family 11, subfamily B, polypeptide 1
Gene ID:
1584

Detection Antibody Information

Catalog Number:
60435-5-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1G2C2
Reactivity:
human
GenBank:
BC096285
Immunogen Catalog Number:
Ag17158

Conjugate:
Unconjugated
Full name:
cytochrome P450, family 11, subfamily B, polypeptide 1
Gene ID:
1584

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50582-3 targets CYP11B1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: CYP11B1 Monoclonal antibody, PBS Only (Capture) 60435-3-PBS (1C8A4). 100 µg. Concentration 1 mg/ml.
Detection antibody: CYP11B1 Monoclonal antibody, PBS Only (Detector) 60435-5-PBS (1G2C2). 100 µg. Concentration 1 mg/ml.
Alternative CYP11B1 matched antibody pairs: MP50582-1, MP50582-2, MP50582-4
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50582-3, CYP11B1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60435-3-PBS. Detection antibody: 60435-5-PBS. Standard: Ag17158. Range: 0.098-50 ng/mL