

For Research Use Only

NSUN4 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50604-1

Capture Antibody Information

Catalog Number:
60448-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1A12E1

Reactivity:
human

Immunogen Catalog Number:
Ag9417

Conjugate:
Unconjugated

Full name:
NOL1/NOP2/Sun domain family, member 4

Gene ID:
387338

Detection Antibody Information

Catalog Number:
60448-2-PBS

Host:
Mouse

Isotype:
IgG2b

Purification Method:
Protein A Magarose purification

Clone ID:
1E9B5

Reactivity:
human

GenBank:
BC014441

Immunogen Catalog Number:
Ag9417

Conjugate:
Unconjugated

Full name:
NOL1/NOP2/Sun domain family, member 4

Gene ID:
387338

Applications

Tested Applications:
Cytometric bead array

Range:
0.195-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50604-1 targets NSUN4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NSUN4 Monoclonal antibody, PBS Only (Capture) 60448-1-PBS (1A12E1). 100 µg. Concentration 1 mg/ml.

Detection antibody: NSUN4 Monoclonal antibody, PBS Only (Detector) 60448-2-PBS (1E9B5). 100 µg. Concentration 1 mg/ml.

Alternative NSUN4 matched antibody pairs: MP50604-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50604-1, NSUN4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60448-1-PBS. Detection antibody: 60448-2-PBS. Standard: Ag9417. Range: 0.195-100 ng/mL.