

For Research Use Only

CNTROB Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50605-2

Capture Antibody Information

Catalog Number:
67061-4-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
1B4E9

Reactivity:
human

Immunogen Catalog Number:
Ag25447

Conjugate:
Unconjugated

Full name:
centrobin, centrosomal BRCA2
interacting protein

Gene ID:
116840

Detection Antibody Information

Catalog Number:
67061-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
3D12A4

Reactivity:
human, mouse, rat, pig

GenBank:
BC021134

Immunogen Catalog Number:
Ag25447

Conjugate:
Unconjugated

Full name:
centrobin, centrosomal BRCA2
interacting protein

Gene ID:
116840

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-100 ng/mL (Cytometric Bead
Array)

Recommended Dilutions:
It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP50605-2 targets CNTROB in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CNTROB Monoclonal antibody, PBS Only (Capture/Detector) 67061-4-PBS (1B4E9). 100 µg. Concentration 1 mg/ml.

Detection antibody: CNTROB Monoclonal antibody, PBS Only (Detector) 67061-1-PBS (3D12A4). 100 µg. Concentration 1 mg/ml.

Alternative CNTROB matched antibody pairs: MP50605-1, MP50605-3, MP50605-4

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50605-2, CNTROB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67061-4-PBS. Detection antibody: 67061-1-PBS. Standard: Ag25447. Range: 1.563-100 ng/mL.