

For Research Use Only

MAOA Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50632-1

Capture Antibody Information

Catalog Number: 68471-2-PBS	Clone ID: 2D7E5	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: monoamine oxidase A
Isotype: IgG1	Immunogen Catalog Number: Ag30028	Gene ID: 4128
Purification Method: Protein G Magarose purification		

Detection Antibody Information

Catalog Number: 68471-3-PBS	Clone ID: 3F1C10	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: monoamine oxidase A
Isotype: IgG1	GenBank: BC008064	Gene ID: 4128
Purification Method: Protein G Magarose purification		
Immunogen Catalog Number: Ag30028		

Applications

Tested Applications: Cytometric bead array	Range: 0.098-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP50632-1 targets MAOA in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MAOA Monoclonal antibody, PBS Only (Capture) 68471-2-PBS (2D7E5). 100 µg. Concentration 1 mg/ml.

Detection antibody: MAOA Monoclonal antibody, PBS Only (Detector) 68471-3-PBS (3F1C10). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50632-1, MAOA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68471-2-PBS. Detection antibody: 68471-3-PBS. Standard: Ag30028. Range: 0.098-100 ng/mL.