

For Research Use Only

# RHOT2 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50638-1

## Capture Antibody Information

**Catalog Number:**  
60464-1-PBS

**Host:**  
Mouse

**Isotype:**  
IgG1

**Purification Method:**  
Protein G Magarose purification

**Clone ID:**  
4B8A11

**Reactivity:**  
human

**Immunogen Catalog Number:**  
Ag1752

**Conjugate:**  
Unconjugated

**Full name:**  
ras homolog gene family, member T2

**Gene ID:**  
89941

## Detection Antibody Information

**Catalog Number:**  
60464-2-PBS

**Host:**  
Mouse

**Isotype:**  
IgG1

**Purification Method:**  
Protein G Magarose purification

**Clone ID:**  
4F7F1

**Reactivity:**  
human

**GenBank:**  
BC014942

**Immunogen Catalog Number:**  
Ag1752

**Conjugate:**  
Unconjugated

**Full name:**  
ras homolog gene family, member T2

**Gene ID:**  
89941

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.098-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50638-1 targets RHOT2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RHOT2 Monoclonal antibody, PBS Only (Capture) 60464-1-PBS (4B8A11). 100 µg. Concentration 1 mg/ml.

Detection antibody: RHOT2 Monoclonal antibody, PBS Only (Capture/Detector) 60464-2-PBS (4F7F1). 100 µg. Concentration 1 mg/ml.

Alternative RHOT2 matched antibody pairs: MP50638-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50638-1, RHOT2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60464-1-PBS. Detection antibody: 60464-2-PBS. Standard: Ag1752. Range: 0.098-100 ng/mL.