

For Research Use Only

# SDHD Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50645-2

## Capture Antibody Information

**Catalog Number:**  
60473-2-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG2a  
**Purification Method:**  
Protein A Magarose purification

**Clone ID:**  
2A12G6  
**Reactivity:**  
human  
**Immunogen Catalog Number:**  
Ag34486

**Conjugate:**  
Unconjugated  
**Full name:**  
succinate dehydrogenase complex, subunit D, integral membrane protein  
**Gene ID:**  
6392

## Detection Antibody Information

**Catalog Number:**  
60473-3-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G Magarose purification

**Clone ID:**  
1G9B8  
**Reactivity:**  
human  
**GenBank:**  
BC012603  
**Immunogen Catalog Number:**  
Ag34486

**Conjugate:**  
Unconjugated  
**Full name:**  
succinate dehydrogenase complex, subunit D, integral membrane protein  
**Gene ID:**  
6392

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.391-200 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50645-2 targets SDHD in immunoassays as a matched antibody pair. Validated in Cytometric bead array.  
Capture antibody: SDHD Monoclonal antibody, PBS Only (Capture/Detector) 60473-2-PBS (2A12G6). 100 µg. Concentration 1 mg/ml.  
Detection antibody: SDHD Monoclonal antibody, PBS Only (Detector) 60473-3-PBS (1G9B8). 100 µg. Concentration 1 mg/ml.  
Alternative SDHD matched antibody pairs: MP50645-1  
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.  
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.  
Antibody use should be optimized for each application and assay.


## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50645-2, SDHD Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60473-2-PBS. Detection antibody: 60473-3-PBS. Standard: Ag34486. Range: 0.391-200 ng/mL.