

For Research Use Only

ARID4B Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50665-4

Capture Antibody Information

Catalog Number:
67384-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
2F9H5

Reactivity:
human

Immunogen Catalog Number:
Ag21485

Conjugate:
Unconjugated

Full name:
AT rich interactive domain 4B (RBP1-like)

Gene ID:
51742

Detection Antibody Information

Catalog Number:
67384-6-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1G4A4

Reactivity:
human

GenBank:
BC130418

Immunogen Catalog Number:
Ag21485

Conjugate:
Unconjugated

Full name:
AT rich interactive domain 4B (RBP1-like)

Gene ID:
51742

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50665-4 targets ARID4B in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ARID4B Monoclonal antibody, PBS Only (Capture) 67384-2-PBS (2F9H5). 100 µg. Concentration 1 mg/ml.

Detection antibody: ARID4B Monoclonal antibody, PBS Only (Detector) 67384-6-PBS (1G4A4). 100 µg. Concentration 1 mg/ml.

Alternative ARID4B matched antibody pairs: MP50665-1, MP50665-2, MP50665-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50665-4, ARID4B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67384-2-PBS. Detection antibody: 67384-6-PBS. Standard: Ag21485. Range: 0.391-100 ng/mL.